



We collect the cells you search for

GILUPI CellCollector®

In vivo Isolation of Circulating Tumor Cells
Detection & Characterization of Circulating Tumor Cells

Imprint


GILUPI GmbH
Am Mühlberg 11
14476 Potsdam
GERMANY

T +49 331 58184782
F +49 331 58184780
E info@gilupi.com

www.gilupi.com
©2015 GILUPI GmbH,
February 2015, V001.

Circulating Tumor Cells

CTCs - a link between primary tumor and metastasis



Source: Klaus Pantel & Ruud H. Brakenhoff, Dissecting the metastatic cascade; Nature Reviews Cancer 4, 448-456 (June 2004)

Cancer related death is usually caused by the spread of invasive tumor cells from the primary tumor and the development of distant metastases. Single tumor cells or cell cluster detach from the primary tumor by down-regulation of cell-cell interactions; they can enter the blood and/or the lymphatic system and will finally leave the bloodstream at an anatomically different location which can lead to the initiation of metastatic growth.¹ The single tumor cells found in the patient's blood are called circulating tumor cells (CTCs).

Liquid biopsy through CTC analysis

Compared to the investigation of the primary tumor or a biopsy taken from a distant metastasis, the investigation of blood is relatively simple, less invasive and can be performed repeatedly. Due to these benefits, CTC investigation can be used in addition to current procedures as a real-time marker for staging, disease progression and therapy responsiveness. Furthermore, the molecular analysis of CTCs can provide highly valuable information for clinicians and researchers as early detection and treatment of the disease and its metastatic spread can reduce cancer mortality dramatically.³

Today, tumor staging is mainly based on tumor size, metastatic lymph node involvement, and evidence of overt distant metastasis obtained by imaging technologies. However, these staging procedures are not sensitive enough to detect early tumor cell spread as a key event in tumor progression. To fill this gap, CTCs have become part of the tumor staging procedure proposed by the American Joint Committee on Cancer (AJCC)² and treatment decisions might be based on these cells in addition to the information obtained from the current procedures.

CTC enumeration and characterization can assist in

- Malignancy detection at an early stage
- Finding treatments tailored to the patient by patient stratification through companion diagnostics
- Treatment monitoring in real time
- Identification of optimal treatments and investigation of drug susceptibility by *ex vivo* culturing of CTCs on an individual level


GILUPI CellCollector®

First *in vivo* platform for circulating tumor cell isolation

CTC detection is challenging due to their extremely low frequency amongst millions of leucocytes and billions of red blood cells. The advanced GILUPI technology aims at the detection and characterization of CTCs combined with the proof of the malignant origin of the collected cells. The tumor cells are captured *in vivo* with high sensitivity and selectivity providing compelling and improved diagnostic value for personalized medicine, early diagnosis as well as treatment and clinical monitoring.

The GILUPI CellCollector® favors “positive enrichment” of CTCs, which is based on the fact that the largest group of cancers is of epithelial origin and cells spread from these tumors express epithelial cell surface markers like EpCAM, the epithelial cell adhesion molecule.⁴ During the *in vivo* application, the functionalized surface comes into direct contact with the circulating blood and rare cells are bound to the device.


In vivo CTC technology - Screening of a larger blood volume




Screening for CTCs in up to 1000ml blood volume

Unique characteristics

GILUPI CellCollector® DC01 (Detektor CANCER01)



GILUPI CellCollector® DC02 (Detektor CANCER02)



Functionalized surface

Tumor cell

Polymer fibers with antibodies


Antibody

The GILUPI CellCollector® is available in two different product versions, Detektor CANCER01 (DC01) recommended for CTC enumeration and Detektor CANCER02 (DC02) recommended for further downstream analysis methods. Both products consist of a stainless steel wire, which has a functionalized surface at one end (gold and hydrogel coating with incorporated anti-EpCAM antibodies). The first generation product DC01 has a smooth, 2 cm long functionalized tip. The next generation product DC02 has a longer, 4 cm functionalized

tip which consists of three thin twisted stainless steel wires for improved haemodynamic characteristics. Both products are CE marked medical devices meeting all requirements of the Medical Devices Directive (MDD) and comply with clinical and diagnostic requirements to be applied in cancer patients. Next-generation products will focus on improved structural characteristics and on the incorporation of different antibodies or antibody combinations for detection of further CTC subtypes (e.g. stem-cell like cells). We also offer customized products.


Workflow

1. *In vivo* CTC isolation in 30 minutes




Product application via standard catheter system

2. Sample preparation



Preparation of isolated cells for subsequent analyses

3. Open for various downstream applications



- Immunofluorescence analysis
- Mutation analysis
- Fluorescence in situ hybridization (FISH)
- *Ex vivo* cell culturing and more

i

Clinical application in solid tumors


Confirmed CTC isolation in

- Lung cancer
- Breast cancer
- Prostate cancer
- Colorectal cancer
- Renal cancer
- Pancreatic cancer
- Neuroendocrine tumors

CTC detection & characterization

- Baseline CTC status
- Before, during and after therapy

Downstream applications




Quantitative detection of CTCs by immunofluorescence analysis

Immunofluorescence microscopy is a powerful technology for determining the presence and subcellular localization of a protein or an antigen in a cell. For this purpose specific fluorescently labeled antibodies are used which bind to the target, and thereby allowing visualization and examination under a fluorescence microscope.

CTCs are identified by positive staining for Keratins and/

or EpCAM and nucleus; negative selection is done via staining for CD45 (marker for blood cells, e.g. leucocytes). Immunofluorescence analysis is successfully employed by GILUPI for quantitative detection of CTCs captured with the GILUPI CellCollector® in various cancer entities like breast, lung, prostate, colorectal and neuroendocrine cancer.



The *in vivo* applied GILUPI CellCollector® DC01 shows a significantly higher CTC detection frequency of 70% in non-metastatic patients and a higher detection frequency of 72% in metastatic patients, compared to the *in vitro* system CELLSEARCH® (CS) (n= 229, pooled data from various cancer entities).


CTCs and their value for targeted therapies

Options to treat cancer patients with individualized, targeted therapy regimens have improved enormously over the last years. Knowledge of the molecular make-up of the cancer cells is an indispensable prerequisite for the effective use of targeted therapies.

Examples for possible downstream applications with the GILUPI CellCollector®:

■ Detection of HER2/neu overexpression in CTCs


In breast cancer, different tumor subtypes vary in terms of biology, histology, therapeutic options, patterns of dissemination and metastasis as well as individual prognostic outcome.⁶ Especially for HER2 positive patients, regular re-testing of remaining cancer cells or CTCs for their HER2 status is important for the adequate choice of targeted therapies (e.g. trastuzumab or lapatinib).



HER2
positive CTC

■ Detection of EML4-ALK translocation in CTCs

Between 2-7% of patients with non-small cell lung cancer (NSCLC) harbor a translocation of the *ALK* gene (European Union).⁷ Lung cancer patients with a proven *EML4-ALK* translocation can benefit from targeted therapies like crizotinib or ceritinib and other *ALK* inhibitors currently under development.⁸ Tumor tissue and CTCs can be tested for the *EML4-ALK* translocation with a companion diagnostic test e.g. Fluorescence in situ hybridization (FISH).



ALK testing by
Fluorescence
in situ hybridization
(FISH)


■ Detection of somatic mutations in CTCs

Identification of mutations in oncogenes (e.g. *EGFR* and *KRAS*) gene can provide information about a patient's responsiveness to targeted therapies (e.g. gefitinib or erlotinib).⁵ Routine testing of CTCs during the treatment can give additional information concerning the usefulness of the applied therapy (e.g. resistance mutation occurrence). CTCs can be characterized on a molecular level using conventional kits. DNA isolated from CTCs can be amplified (whole genome amplification) and tested for therapy relevant mutations using e.g. mutation specific quantitative real-time PCR. In this example, a *KRAS* Mutation is detected in a lung cancer patient.

Further applications

- Ex vivo cell culturing
- Next generation sequencing (NGS)
- Gene expression analysis
- Protein analysis

KRAS Mutation Analysis



GILUPI CellCollector®



CTCs isolated with the GILUPI CellCollector® technology can be enumerated and characterized to assist in

- Malignancy detection at an early stage
- Finding treatments tailored to the patient by patient stratification through companion diagnostics
- Treatment monitoring in real time
- Identification of optimal treatments and investigation of drug susceptibility by *ex vivo* culturing of CTCs on an individual level

i

Applications

- CTC enumeration
- Immunofluorescence analysis of proteins of interest (e.g. HER2, ALK, AR-V7)
- Mutation analysis (e.g. *EGFR*, *KRAS*)
- Fluorescence in situ hybridization (FISH)
- Next generation sequencing (NGS)

Further applications are currently established (e.g. mRNA analysis)

GILUPI offers



■ Diagnostic Services

Quantitative detection of CTCs (immunofluorescence analysis) in an in-house diagnostics lab



■ Customized products

Expertise in incorporation of antibodies or antibody combinations on customized GILUPI CellCollector®

Detection of further CTC subtypes (by other antibodies than EpCAM)

References

- 1 Pantel, K., Brakenhoff, R.H. (2004). Dissecting the metastatic cascade. *Nature Reviews Cancer* 4, 448-456.
- 2 Edge, S.B. and Compton, C.C. (2010). The American Joint Committee on Cancer: the 7th edition of the AJCC cancer staging manual and the future of TNM. *Annals of surgical oncology* 17, 1471-1474.
- 3 Plaks, V., Koopman, C.D., and Werb, Z. (2013). Cancer. Circulating tumor cells. *Science* 341, 1186-1188.
- 4 Saucedo-Zeni, N., Mewes, S., Niestroj, R., Gasiorowski, L., Murawa, D., Nowaczyk, P., Tomasi, T., Weber, E., Dworacki, G., Morgenthaler, N.G., Jansen, H., Propping, C., Sterzynska, K., Dyszkiewicz, W., Zabel, M., Kiechle, M., Reuning, U., Schmitt, M., and Lucke, K. (2012). A novel method for the in vivo isolation of circulating tumor cells from peripheral blood of cancer patients using a functionalized and structured medical wire. *Int. J. Oncol.* 41, 1241-1250.
- 5 Suda, K., Tomizawa, K., and Mitsudomi, T. (2010). Biological and clinical significance of *KRAS* mutations in lung cancer: an oncogenic driver that contrasts with *EGFR* mutation. *Cancer Metastasis Rev.* 29, 49-60.
- 6 Perou, C.M., Sorlie, T., Eisen, M.B., van de Rijn, M., Jeffrey, S.S., Rees, C.A., Pollack, J.R., Ross, D.T., Johnsen, H., Akslen, L.A., Fluge, O., Pergamenschikov, A., Williams, C., Zhu, S.X., Lonning, P.E., Borresen-Dale, A.L., Brown, P.O., and Botstein, D. (2000). Molecular portraits of human breast tumours. *Nature* 406, 747-752.
- 7 Scagliotti, G., Stahel, R.A., Rosell, R., Thatcher, N., and Soria, J.C. (2012). *ALK* translocation and crizotinib in non-small cell lung cancer: an evolving paradigm in oncology drug development. *Eur. J. Cancer* 48, 961-973.
- 8 Peters, S., Taron, M., Bubendorf, L., Blackhall, F., and Stahel, R. (2013). Treatment and detection of *ALK*-rearranged NSCLC. *Lung Cancer* 81, 145-154.



GILUPI CellCollector® is available for distribution and sale in countries accepting CE mark.
GILUPI CellCollector® is available for research purposes (*in vitro* only) worldwide.
For further information on intended use, warnings and limitations please refer to the GILUPI CellCollector®
Instructions for Use, our homepage or contact us directly.

Visit gilupi.com

GILUPI GmbH

Am Mühlenberg 11
14476 Potsdam
GERMANY

T +49 331 58184782
F +49 331 58184780
E info@gilupi.com